Technical support: order@acebiolab.com

Phone: 886-3-2870051

Ver.1 Date: 20220104

AceColor Prestained Protein Marker (25-400 kDa)

Cat# A1058 - 2 x250 μ l Storage at -20°C

INFORMATION

ACE Biolabs AceColor Prestained Protein Ladder is a mixture of 7 blue-, 2 orange-, and a green-stained proteins (25 to 400 kDa) for use as size standards in protein electrophoresis (SDS-PAGE) and Western blotting. The protein ladder is supplied in a ready-to-use format for direct loading onto gels; no need to heat, reduce, or add sample buffer prior to use.

APPLICATIONS

- Monitoring protein migration during SDS-PAGE.
- Monitoring protein transfer efficiency on membranes after Western blot transferring.
- Sizing estimation of proteins on SDS-PAGE and Western blots.

PROTOCOL

- For SDS-PAGE: Loading 3 μl or 5 μl per well for clear monitoring electrophoresis on a 15-well or 10-well mini-gel. Apply more for a thicker (> 1.5 mm) or larger gel.
- For Western blotting: Loading 3~5 μl per well for monitoring electrophoresis and general Western blot transferring efficiency. Apply more for a thicker (> 1.5 mm) or larger gel.

Note: Do NOT heat, dilute, or add reducing agents before loading.

LADDER

Gel type Gel concentration Running buffer		Tris-Glycine			Tris-Acetate		Bis-Tris	
		6%	7%	B4-20%	6%	T3-8%	T4-12%	T4-12%
		Tris-Glycine			Tris-Acetate		MES	MOPS
		Apparent Molecular Weights, kDa						
	10		400	0.000	350		350	350
	20	400 300 250	400 300 250 200 160	350 300 250 200 160	350 300 250 200		350 300 200 200 160 130 100	350 300 200 160 130
	30	250 200 160	130	= 200 160 130	— 160 — 130	= 350 300 250	20000	- 130 - 100
gel	40	130	100	 100	100	200 160	 65	65
% length of gel	50	100	— 72	 70		- 130	— 45	— 65
engt	60			— 45	 65	— 100		— 45
%1	70	 72	— 45	— 25	— 45	— 65	— 25	
	80		45			— 45		— 25
	90	— 45 — 25						
	100		— 25		— 25	— 25		

PRODUCT USE LIMITATION

These products are intended for research use only.

