

Ponceau S solution

Cat# A1052

Store at 2 to 8°C

INTRODUCTION

It is commonly used for the detection of protein on PVDF membrane, nitrocellulose membrane and cellulose acetate membrane; immerse the PVDF membrane, nitrocellulose membrane or cellulose acetate membrane with protein in the Ponceau staining solution and shake for 5-10 Minutes or more, after the protein bands appear on the membrane, take out the blotting membrane and rinse with distilled water, PBS or other appropriate solution for 2-3 times, each time for 3-5 minutes, until the red background disappears, the subsequent follow-up WB detection.

PRODUCT COMPONENTS

Acetic acid	5%
Ponceau S	0.1%

PRODUCT USE LIMITATION

These products are intended for research use only.