

# 2X ACE LAMP Master Mix

Cat# EP1407 – 1ml | EP1408 – 5\*1 ml | EP1409 – 15\*1 ml

Storage: All components should be stored at -20°C.

## INTRODUCTION

2× ACE LAMP Master Mix is a blend of Taq DNA Polymerase and a DNA proofreading polymerase with 3' to 5' exonuclease activity. Its fidelity was 6-fold higher than conventional Taq DNA Polymerase. Used with the optimized buffer system, 2× ACE LAMP Master Mix is applicable to long PCR products, up to 21 kb. This Master Mix is also able to amplify long fragments accurately from templates of different sources or different length.

2× ACE LAMP Master Mix contains Vazyme LAMP DNA Polymerase, dNTP, and optimized buffer. The reaction can be started by adding only primers and template, which simplifies the operation, improves through-put, and enhances result reproducibility. The protective agents included guarantees the stability of the activity of this Master Mix. The PCR product, containing dA at 3'-end, can be cloned into T-vector, and is suitable for One Step Express cloning kit.

## CONTENTS

Component	EP1407	EP1408	EP1409
2X ACE LAMP Master Mix	1ml	5* 1 ml	15* 1 ml

## PROTOCOL

### 1. General reaction mixture for PCR:

ddH <sub>2</sub> O	to 50 µl
2× ACE LAMP Master Mix	25 µl
Template DNA*	Optional
Primer 1 (10 µM)	2 µl
Primer 2 (10 µM)	2 µl

\*The recommended amount of DNA template for a 50 µl reaction system is as follows:

Human Genomic DNA	10 - 200 ng
Bacterial Genomic DNA	1 - 100 ng
λ DNA	0.1 - 10 ng
Plasmid DNA	0.1 - 10 ng

## 2. Thermocycling conditions:

---

94°C	5 min (Pre-denaturation)	
94°C	30 sec	} 30 - 35 cycles
55°C *	30 sec	
72°C	30 sec / kb	
72°C	7 min (Final extension)	

---

\*The optimal annealing temperature should be 1-2°C lower than the  $T_m$  of the primers used.

---

94°C	1-3 min (Pre-denaturation)	
94°C	10 sec	} 30 - 35 cycles
68°C *	30-60sec/kb	
68°C	7 min (Final extension)	

---

\* For amplification of a DNA fragment > 5 kb, it is recommended to use long primers which  $T_m$  between 68°C and 70°C. The temperature for both annealing and extension should be 68°C, which can significantly improve the amplification specificity. Extending extension time could increase the amplification yield.

## PRIMERS DESIGNING NOTES

1. Choose C or G as the last base of the 3'-end of the primer;
2. Avoid continuous mismatching at the last 8 bases of the 3'-end of the primer;
3. Avoid hairpin structure at the 3'-end of the primer;
4.  $T_m$  of the primers should be within the range of 55°C - 65°C;
5. Additional sequence should not be included when calculating  $T_m$  of the primers;
6. GC content of the primers should be within the range of 40% - 60%;
7.  $T_m$  and GC content of forward and reverse primers should be as similar as possible.

## PRODUCT USE LIMITATION

These products are intended for research use only.