

p38 (phospho Thr180/Y182) Polyclonal Antibody

 Cat# A20311PI
 store at -20°C

INFORMATION

Size	50 ul
Immunogen	The ant
Reactivity	Human,Mouse,Rat,Guineapig
Applications	IF,WB,IHC-p,ELISA
Immunogen	The antiserum was produced against synthesized peptide derived from human p38 MAPK around the phosphorylation site of Thr179 and Tyr181. AA range:151-200
Specificity	Phospho-p38 (T180/Y182) Polyclonal Antibody detects endogenous levels of p38 protein only when phosphorylated at T180/Y182.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Dilution	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Source	Rabbit
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage Stability	-20°C/1 year
Concentration	1 mg/ml
Molecular Weight	38KD







