


# Hemin

Cat# C45210B  
store at 2-8°C

## INFORMATION

Size	1g
Description	Iron-containing compound of porcine origin that can induce hydrogen peroxide-dependent lipid peroxidation.
Synonym(s)	Hematin chloride ; Hemin ; Hematine hydrochloride ; Chloroproporphyrin IX iron(III) ; 1,3,5,8-Tetramethyl-2,4-divinylporphine-6,7-dipropionic acid ferrichloride ; Chlorohemin ; Chloroferriporphyrin ; Ferriheme chloride ; Ferriporphyr in chloride
Biological source	Porcine
CAS#	16009-13-5
M.W.	651.96 (C <sub>34</sub> H <sub>32</sub> C <sub>1</sub> N <sub>4</sub> FeO <sub>4</sub> )
Purity	≥98.0%
Storage temp	2-8°C
Solubility	1.4 M NaOH: 25mg/mL



## PRODUCT USE LIMITATION