

Rabbit Anti-Rat IgG (H+L), HRP

Cat# A120111
store at -20°C

INFORMATION

Size	1 ml
Host/Isotope	Rabbit
Reactivity	rabbit
Applications	ELISA,WB,IHC
Category	Antigen antibody
Antibody concentration	1-1.2mg/ml
Purification method	Antigen affinity purification
Buffer	0.02M PBS, 50% glycerol, pH7.4
Stabilizer	0.2%BSA
Preservative	0.02% thimerosal
Use concentration	ELISA: 1/4000-1/20000; WB: 1/4000-1/20000; IHC: 1/500-1/1000; Note: Other combinations can be used according to actual Proper proportion
Storage Stability	Transport at 4°C, store at -20°C away from light
Precautions	<ol style="list-style-type: none">1. In order to obtain the best results of the experiment, please adjust the best working dilution according to the experiment.2. When using, it should be diluted according to the required concentration, and ready to use.3. If it is seldom used at a time, it can be used after appropriate aliquots, which can effectively prevent reagent contamination.4. Although this reagent is in you can still use it for a long time after the validity period, but it is strongly recommended that you use it correctly within the validity period specified by the tube.5. This product is only used for scientific research and cannot be used for human experiments or human treatment..