




Human IP-10 Polyclonal Antibody

Cat# A340006

store at Store at -20°C. Avoid freeze / thaw cycles.

INFORMATION

Size	100 ug
Applications	Antibody
Host	Rabbit
Note	Centrifuge before opening to ensure complete recovery of vial contents.
Synonyms	CXCL10,C7,IFI10,INP10,IP-10,SCYB10,crg-2,gIP-10,mob-1
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000
Storage Stability	Store at -20°C. Avoid freeze / thaw cycles.
Background	This antimicrobial gene encodes a chemokine of the CXC subfamily and ligand for the receptor CXCR3. Binding of this protein to CXCR3 results in pleiotropic effects, including stimulation of monocytes, natural killer and T-cell migration, and modulation of adhesion molecule expression.



PRODUCT USE LIMITATION

These products are intended for research use only.