




Human CXCL6 Antibody

Cat# A340008

store at Store at -20°C. Avoid freeze / thaw cycles.

INFORMATION

Size	100 ug
Applications	Antibody
Host	Rabbit
Note	Centrifuge before opening to ensure complete recovery of vial contents.
Synonyms	CKA 3,CKA-3,CKA3,CXCL 6,CXCL6,CXCL6,GCP2,GCP-2,Granulocyte chemotactic protein 2
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4
Purify	Antigen affinity purification
Dilution	IHC 1:40-1:200, ELISA 1:5000-1:10000
Background	Chemotactic for neutrophil granulocytes. Signals through binding and activation of its receptors (CXCR1 and CXCR2). In addition to its chemotactic and angiogenic properties, it has strong antibacterial activity against Gram-positive and Gram-negative bacteria (90-fold-higher when compared to CXCL5 and CXCL7).



Immunohistochemistry of paraffinembedded
Human thyroid cancer tissue using CXCL6
Polyclonal Antibody at dilution of 1:45(×200)

PRODUCT USE LIMITATION

These products are intended for research use only.