




(-)-Galocatechin

Cat# C016092S

store at at 2-8°C.

INFORMATION

Size	20mg
Applications	Chemicals>Natural Products
M.F	C15H14O7
Physical Description	White powder
Purity	> 98%
M.W	306.27
Synonyms	(2S)-3,4-Dihydro-2-(3,4,5-trihydroxyphenyl)-2H-1-benzopyran-3,5,7-triol.
Storage Stability	2-8°C, Protected from air and light, refrigerate or freeze.
Solvent	Pyridine, Alcohol, Methanol, H2O, etc.
Source	The wood of Acacia catechu (L.F.) Willd.
HPLC Method	Mobile phase: Methanol - 0.05%Trifluoroacetic acid in water, gradient elution; Flow rate: 1.0 ml/min;Column temperature: Room Temperature;The wave length of determination: 210 nm.
Intended Use	1. Reference standards; 2. Pharmacological research; 3. Synthetic precursor compounds; 4. Intermediates & Fine Chemicals; 5. Others.
Biological Activity or Inhibitors	Galocatechins have antioxidant potential.



PRODUCT USE LIMITATION