


KIT Antibody

Cat# A0118P

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

INFORMATION

Product Name	KIT Antibody	Western Blot analysis of 453 cells using c-Kit Polyclonal Antibody 
Cat. No.	A0118P	
Uniprot No.	P10721	
Product type	Primary antibodies	
Species Reactivity	Human , Mouse	
Host	Rabbit	
Clonality	Polyclonal	
Immunogen	Synthesized peptide derived from human c-Kit around the non-phosphorylation site of Y936.	
Tested applications	WB , IHC , ELISA	
Application	WB : 1 : 500 - 2000 IHC : 1 : 100 – 300 ELISA : 1 : 5000	
Conjugation	Non-conjugated	
Purification method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.	
Storage buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.	
Storage instruction	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.	
Alias	KIT; SCFR; Mast/stem cell growth factor receptor Kit; SCFR; Piebald trait protein; PBT; Proto-oncogene c-Kit; Tyrosine-protein kinase Kit; p145 c-kit; v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog; CD antigen CD117	

PRODUCT USE LIMITATION

These products are intended for research use only.