

Goat Anti-Mouse IgG (H+L)

Cat# A1010S – 1mg / A1010S – 10mg

Storage at 2-8°C

INFORMATION

Product Name	Goat Anti-Mouse IgG (H+L)
Cat. No.	A1010S
Product type	Secondary antibodies
Species Reactivity	Mouse
Immunogen	Mouse IgG, whole molecule
Host	Goat
Clonality	Polyclonal
Target	Mouse IgG
Tested applications	WB, ICC/IF, ELISA
Application	WB : 1:5000-10000 ELISA : 1:5000-10000
Conjugation	HRP
Purification	Purified polyclonal antibody from Goat serum with protein A column
Storage buffer	Phosphate Buffered Saline, pH7.4
Storage instruction	2-8°C

PRODUCT USE LIMITATION

These products are intended for research use only.