

# Goat Anti-Rabbit IgG (H+L)

Cat# A1011S – 1mg / A1011S – 10mg

Storage at 2-8°C

## INFORMATION

<b>Product Name</b>	Goat Anti-Rabbit IgG (H+L)
<b>Cat. No.</b>	A1011S
<b>Product type</b>	Secondary antibodies
<b>Species Reactivity</b>	Rabbit
<b>Immunogen</b>	Rabbit IgG, whole molecule
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Target</b>	Rabbit IgG
<b>Tested applications</b>	WB, ICC/IF, ELISA
<b>Application</b>	WB : 1:5000-10000 ELISA : 1:5000-10000
<b>Conjugation</b>	HRP
<b>Purification</b>	Purified polyclonal antibody from Goat serum with protein A column
<b>Storage buffer</b>	Phosphate Buffered Saline, pH7.4
<b>Storage instruction</b>	2-8°C

## PRODUCT USE LIMITATION

These products are intended for research use only.