

Mouse anti-Human IgG (Fab)

Cat# A120105– 1 mg

Storage at -20 °C

INFORMATION

Product Name	Mouse anti-Human IgG (Fab)
Cat. No.	A120105
Product type	Secondary antibodies
Species Reactivity	Human
Host	Mouse
Target	Human IgG
Size	1 mg
Tested applications	WB, ELISA
Application	Optimal concentrations should be determined by the end user
Conjugation	None
Purification	Affinity purified with protein A or protein G of ascites
Dilute buffer	0.01M PBS, pH7.4
Storage instruction	Upon delivery aliquot and store at -20 °C or -80 °C . Avoid repeated freeze/thaw cycles.

PRODUCT USE LIMITATION

These products are intended for research use only.