

# Goat anti-Mouse IgG (H+L) for Control line

Cat# A120108– 50 mg

Storage at -20 °C

## INFORMATION

<b>Product Name</b>	Goat anti-Mouse IgG (H+L) for Control line
<b>Cat. No.</b>	A120108
<b>Product type</b>	Secondary antibodies
<b>Species Reactivity</b>	Mouse
<b>Host</b>	Goat
<b>Target</b>	Mouse IgG
<b>Size</b>	50 mg
<b>Tested applications</b>	WB, ELISA
<b>Application</b>	Optimal concentrations should be determined by the end user
<b>Conjugation</b>	None
<b>Purification</b>	Affinity purified
<b>Dilute buffer</b>	0.01M PBS, pH7.4
<b>Storage instruction</b>	Upon delivery aliquot and store at -20 °C or -80 °C . Avoid repeated freeze/thaw cycles.

## PRODUCT USE LIMITATION

These products are intended for research use only.