

Goat anti-Rabbit IgG (H+L) for Control line

Cat# A120109– 50 mg

Storage at -20 °C

INFORMATION

Product Name	Goat anti-Rabbit IgG (H+L) for Control line
Cat. No.	A120109
Product type	Secondary antibodies
Species Reactivity	Rabbit
Host	Goat
Target	Rabbit IgG
Size	50 mg
Tested applications	WB, ELISA
Application	Optimal concentrations should be determined by the end user
Conjugation	None
Purification	Affinity purified
Dilute buffer	0.01M PBS, pH7.4
Storage instruction	Upon delivery aliquot and store at -20 °C or -80 °C . Avoid repeated freeze/thaw cycles.

PRODUCT USE LIMITATION

These products are intended for research use only.