

# Rabbit Anti-Bovine IgG (H+L)

Cat# A28637PI – 100 µl

Storage at -20 °C

## INFORMATION

<b>Product Name</b>	Rabbit Anti-Bovine IgG (H+L)
<b>Cat. No.</b>	A28637PI
<b>Product type</b>	Secondary antibodies
<b>Species Reactivity</b>	Bovine
<b>Immunogen</b>	Bovine IgG, whole molecule
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Target</b>	Bovine IgG
<b>Concentration</b>	4-5 mg/ml
<b>Tested applications</b>	WB, ICC/IF, IHC-p, ELISA
<b>Application</b>	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000, IF 1:5000-10000
<b>Conjugation</b>	None
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Storage buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage instruction</b>	Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

## PRODUCT USE LIMITATION

These products are intended for research use only.