


Anti-Rabies virus Glycoprotein Antibody

Cat# A05003PC

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Anti-Rabies virus Glycoprotein Antibody	
Cat. No.	A05003PC	
Size	50 µg, 100 µg	
Uniprot	P03524	
Product type	Primary antibodies	
Species Reactivity	Rabies virus	
Immunogen	Recombinant Rabies virus Glycoprotein protein (20-459AA)	
Host	Rabbit	
Clonality	Polyclonal	
Tested applications	ELISA, WB	
Application	WB:1:500-1:5000	
Conjugation	Non-conjugated	
Purification Method	Protein G purified, >95%	
Storage buffer	50% Glycerol, 0.01M PBS, PH 7.4	
Alias	Glycoprotein	
Image		Western blot All lanes: Rabies virus Glycoprotein antibody at 1: 500 Lane 1: Rabies virus antigen 1: 2 Lane 2: Rabies virus antigen 1: 10 Lane 3: Rabies virus antigen 1: 20 Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution

PRODUCT USE LIMITATION

These products are intended for research use only.