

# Polyclonal Goat Anti-Mouse IgG2a, HRP

Cat# A28751PI

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	Polyclonal Goat Anti-Mouse IgG2a, HRP
<b>Cat. No.</b>	A28751PI
<b>Size</b>	100 ul
<b>Product type</b>	Secondary antibody
<b>Reactivity</b>	MouseIgG2a
<b>Host</b>	Goat
<b>Concentration</b>	4-5mg/mL
<b>Tested applications</b>	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
<b>Application</b>	WB,IHC-p,IF,ELISA
<b>Conjugation</b>	HRP
<b>Purification Method</b>	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
<b>Storage buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## PRODUCT USE LIMITATION

These products are intended for research use only.