


Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) mouse mAb

Cat# A33261PI

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

INFORMATION

Product Name	Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) mouse mAb
Cat. No.	A33261PI
Size	50ug, 100ug
Product type	Primary antibody
Species Reactivity	Hu, Ms
Immunogen	Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) synthetic peptide conjugated to KLH.
Host	Mouse
Clonality	Monoclonal
Application	WB
Concentration	1 mg/ml
Conjugation	Non-conjugated
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Alias	ERK 1;ERK 2;ERK-2;ERK1;erk1/2;ERK2;ERT1;ERT2;Extracellular signal regulated kinase 1;Extracellular signal regulated kinase 1;Extracellular signal regulated kinase 2;Extracellular signal regulated kina



PRODUCT USE LIMITATION

These products are intended for research use only.