

Integrin Antibody Sampler Kit

Cat# AK0178

Upon receipt, store at -20°C. Avoid freeze/thaw cycles.

PRODUCT DESCRIPTION

Integrins are α/β heterodimeric cell surface receptors that play a pivotal role in cell adhesion and migration, as well as in growth and survival. The integrin family contains at least 18 α and 8 β subunits that form 24 known integrins with distinct tissue distribution and overlapping ligand specificities. Integrins not only transmit signals to cells in response to the extracellular environment (outside-in signaling), but also sense intracellular cues to alter their interaction with the extracellular environment (inside-out signaling).

PRODUCT INCLUDES

Cat No.	Product name	Quantity	Applications	Reactivity	Host
A340723	ITGA5 Polyclonal Antibody	20 μ L	WB, ELISA	Human, Mouse, Rat	Rabbit
A340563	ITGAV Polyclonal Antibody	20 μ L	WB, IHC, ELISA	Human, Mouse, Rat	Rabbit
A340716	ITGA4 Polyclonal Antibody	20 μ L	WB, IHC, ELISA	Human, Mouse, Rat	Rabbit
A340564	ITGB1 Polyclonal Antibody	20 μ L	WB, ELISA	Human, Mouse, Rat	Rabbit
A340565	CD61 Polyclonal Antibody	20 μ L	WB, IHC, ELISA	Human, Mouse, Rat	Rabbit
A340566	ITGB5 Polyclonal Antibody	20 μ L	WB, ELISA	Human, Mouse, Rat, Monkey	Rabbit
A1013s	Goat Anti-Rabbit IgG (H+L) (peroxidase/HRP conjugated)	120 μ L	WB, ELISA	Rabbit	Goat

PRODUCT USE LIMITATION

These products are intended for research use only.