

Jak Isoform Sampler Kit

Cat# AK0181

Upon receipt, store at -20°C. Avoid freeze/thaw cycles.

PRODUCT DESCRIPTION

Members of the Janus family of tyrosine kinases (Jak1, Jak2, Jak3, and Tyk2) are activated by ligands binding to a number of associated cytokine receptors. Upon cytokine receptor activation, Jak proteins become autophosphorylated and phosphorylate their associated receptors to provide multiple binding sites for signaling proteins. These associated signaling proteins, such as Stats, Shc, insulin receptor substrates, and focal adhesion kinase (FAK), typically contain SH2 or other phospho-tyrosine-binding domains.

PRODUCT INCLUDES

Cat No.	Product name	Quantity	Applications	Reactivity	Host
A340575	JAK1 Polyclonal Antibody	20µL	WB, IHC, ELISA	Human, Mouse, Rat	Rabbit
A340576	JAK2 Polyclonal Antibody	20µL	WB, IHC, ELISA	Human, Mouse, Rat	Rabbit
A340577	JAK3 Polyclonal Antibody	20µL	WB, IHC, IF, ELISA	Human, Mouse, Rat	Rabbit
A1013s	Goat Anti-Rabbit IgG (H+L) (peroxidase/HRP conjugated)	120µL	WB, ELISA	Rabbit	Goat

PRODUCT USE LIMITATION

These products are intended for research use only.