

# Phospho-EGF Receptor Pathway Antibody Sampler Kit

Cat# AK0224

Upon receipt, store at -20°C. Avoid freeze/thaw cycles.

## PRODUCT DESCRIPTION

The epidermal growth factor (EGF) receptor is a transmembrane tyrosine kinase that belongs to the HER/ErbB protein family. Ligand binding results in receptor dimerization, autophosphorylation, activation of downstream signaling, internalization, and lysosomal degradation. Phosphorylation of EGF receptor (EGFR) at Tyr845 in the kinase domain is implicated in stabilizing the activation loop, maintaining the active state enzyme, and providing a binding surface for substrate proteins. c-Src is involved in phosphorylation of EGFR at Tyr845. The SH2 domain of PLC  $\gamma$  binds at phospho-Tyr992, resulting in activation of PLC  $\gamma$ -mediated downstream signaling. Phosphorylation of EGFR at Tyr1045 creates a major docking site for the adaptor protein c-Cbl, leading to receptor ubiquitination and degradation following EGFR activation. The GRB2 adaptor protein binds activated EGFR at phospho-Tyr1068. A pair of phosphorylated EGFR residues (Tyr1148 and Tyr1173) provide a docking site for the Shc scaffold protein, with both sites involved in MAP kinase signaling activation. Phosphorylation of EGFR at specific serine and threonine residues attenuates EGFR kinase activity. EGFR carboxy-terminal residues Ser1046 and Ser1047 are phosphorylated by CaM kinase II; mutation of either of these serines results in upregulated EGFR tyrosine autophosphorylation.

## PRODUCT INCLUDES

Cat No.	Product name	Quantity	Applications	Reactivity	Host
<b>A340305</b>	Phospho-EGFR (Tyr1069) Polyclonal Antibody	20 $\mu$ L	WB, ELISA	Human, Mouse, Rat, Monkey	Rabbit
<b>A340365</b>	Phospho-AKT1 (Ser473) Polyclonal Antibody	20 $\mu$ L	WB, IHC, IF, ELISA	Human, Mouse, Rat	Rabbit
<b>A340228</b>	Phospho-GAB1 (Tyr627) Polyclonal Antibody	20 $\mu$ L	WB, ELISA	Human, Mouse, Monkey	Rabbit
<b>A340370</b>	Phospho-ERK 1/2 (Thr202/Tyr204) Polyclonal Antibody	20 $\mu$ L	WB, IHC, ELISA	Human, Mouse, Rat	Rabbit
<b>A340324</b>	Phospho-SHC1 (Tyr349) Polyclonal Antibody	20 $\mu$ L	WB, IHC, ELISA	Human, Mouse, Rat	Rabbit
<b>A340279</b>	Phospho-PLCG1 (Tyr783) Polyclonal Antibody	20 $\mu$ L	WB, IHC, ELISA	Human, Mouse, Rat, Monkey	Rabbit

<b>A1013s</b>	Goat Anti-Rabbit IgG (H+L) (peroxidase/HRP conjugated)	120μL	WB, ELISA	Rabbit	Goat
---------------	---	-------	-----------	--------	------

## **PRODUCT USE LIMITATION**

These products are intended for research use only.