

## Kinetin (6KT)

Cat# C3311

Store at Room temperature

### INTRODUCTION

**CAS:** 525-79-1

**M.W:** 215.21

**Formula:** C<sub>10</sub>H<sub>9</sub>N<sub>5</sub>O

**Melting Point :** 265 ~ 266°C


**Purity:** ≥99.0%(HPLC)

**Appearance:** White crystalline powder

**Solubility:** Easily soluble in acetone and ether, slightly soluble in chloroform, insoluble in water

**Synonym :** 6-Furfuryl Aminopurine ; Kinetin ; 6-KT

### CHEMICAL STRUCTURE



### PRODUCT USE LIMITATION

These products are intended for research use only.