

Trisodium citrate dihydrate


Cat# C3333

Store at room temperature

INTRODUCTION

Grade: AR**CAS:** 6132-04-3**M.W:** 294.1**Formula:** $C_6H_5Na_3O_7 \cdot 2H_2O$ **Purity:** 99%**Solubility:** Soluble in water (29.4 g/l) at 20 °C. Insoluble in ethanol**Density:** 1.76**Melting Point:** 300°C**Synonym:** Citric acid trisodium salt dihydrate

CHEMICAL STRUCTURE



PRODUCT USE LIMITATION

These products are intended for research use only.