

Pepstatin A

Cat# C4086

Store at 2-8 °C


INTRODUCTION

CAS: 26305-03-3**M.W:** 685.90**Formula:** C₃₄H₆₃N₅O₉**Purity:** ≥75% (HPLC)**Melting Point:** 220~240°C**Specific rotation:** -85 ~-95° (0.5%, MeOH)**Appearance:** White powder**Solubility:** 9 : 1 = methanol : acetic acid : 1 mg/ml ; Ethanol : 1-2 mg/ml with heat up to 60°C ; DMSO : 5 mg/ml ;

Insoluble in benzene, chloroform, water, 1 M NaOH and ether.

Synonym : Pepstatin

CHEMICAL STRUCTURE



PRODUCT USE LIMITATION

These products are intended for research use only.