

pLV-ACE2-3xFLAG-IRES-puro

Cat# PV010

Store at -20°C


INTRODUCTION

Vector : Lenti**Prokaryote resistance:** AMP⁺**Expression type:** OE**Plasmid tag:** Puromycin**Cloning site:** EcoR I, BamH I**Tag start sequence:** IRES**Promoter:** CMV**Primer:**

CMV-Forward: 5'-CGCAAATGGGCGGTAGGCGTG-3'

G494-S-IRES-N-3 : 5'-GGTTCCGGGCCCTCACATTG-3'

VECTOR INFORMATION



PRODUCT USE LIMITATION

These products are intended for research use only.