

10X PBS Buffer, sterile

Cat# C1015

Store at room temperature, maintain sterility

INTRODUCTION

PBS (phosphate buffered saline) is the most widely used balanced salt solution in biological research, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting, and preparing reagents. It is diluted to a 1X working concentration

Product information

| | |
|----------------------|--------------------------|
| Form | Clear, colorless liquid |
| Product size | 1L |
| pH range | 7.3-7.5 |
| Inorganic Salts | No Calcium, No Magnesium |
| Phenol Red Indicator | No Phenol Red |

Product components

| | |
|----------------------------------|-------|
| NaCl | 37mM |
| KCl | 2.7mM |
| Na ₂ HPO ₄ | 8mM |
| KH ₂ PO ₄ | 2mM |

PRODUCT USE LIMITATION

These products are intended for research use only.