

Polyethyleneimine

Cat# C2020

Storage at Room Temperature

INFORMATION

CAS: 106899-94-9**M.W.:** 600**Resin cont (. wt%):** >98**Viscosity (mPa · s-25°C):** 500-2500**pH (5%aq.):** 10-12**Appearance:** Colorless or light yellow liquid**Specific gravity:** 1.03**Amine value (mmol / g, solid):** 20**Freezing point:** <-20 °C**Decomposition Temperature:** 270 °C**Flash point:** 263 °C**Solubility:** Soluble in water, alcohol. Partly soluble in ethyl acetate, THF, toluene.**Aminerate:**

Primary - 35%

Secondary - 35%

Tertiary - 30%

PRODUCT USE LIMITATION

These products are intended for research use only.