

# Agar

Cat# C3065

Storage at room temperature

## DESCRIPTION


For the solidification of culture medium, immobilized enzyme carriers and gel electrophoresis.

## INFORMATION

**Chemical Formula** :  $(C_{12}H_{18}O_9)_N$ **CAS No.** : 9002-18-0**Molecular Weight** : 336.33496**Solubility** : Soluble in boiling water

**Alias** : agar-epi; macconkey agar; macconkey agar cs; mac conkey agar no 1; macconkey agar no 1, vegitone; macconkey agar with mug; macconkey agar without crystal violet; macconkey bromocresol purple broth; macconkey broth purple; mac conkey mug agar; vegitone macconkey agar no 1; agaragarflake; agar-agargum; agaropectin, mixtwithagarose; agarose, mixtwithagaropectin; bengalgelatin; bengalisinglass; ceylon; ceylonisinglass

## CHEMICAL STRUCTURES



## PRODUCT USE LIMITATION

These products are intended for research use only.