

M17 Broth

Cat# C3082 – 250g

Storage at Room Temp and dry environment

APPLICATION

For culture *lactic streptococci* and their bacteriophages.

FORMULA

Component	g/L
Tryptone	5g
Soy Peptone	5g
Yeast extract	2.5g
Meat extract	5g
Ascorbic acid	0.5g
Magnesium sulfate	0.25g
Di-sodium-glycerophosphate	19g
Lactose	5g
pH 6.9 ± 0.2 @ 25°C	

PROTOCOL

1. Weight 42.3 g M17 Broth and mix with 1000 ml ddH₂O thoroughly.
2. Sterilize by autoclaving at 121°C for 15 mins.
3. Store at 2-8°C. If there is any physical changes, discard the medium.

PRODUCT USE LIMITATION

These products are intended for research use only.