

Normal Goat Serum

Cat# P0004GtS

Storage at -20°C for 5 years

INFORMATION

Product Name	Normal Goat Serum
Product Type	Serum
Cat NO.	P0004GtS
Size	10 / 50 ml
Description	The Normal Goat Serum is collected from healthy goat by aseptic collection and ensures free of bacteria, fungi, mycoplasma and virus contamination by passing three 0.1µm filters. All operations comply with the GMP (Good Manufacturing Practices for Drug) Production Standards.
Application	Cell culture, blocking step of immunofluorescence staining and immunohistochemistry et al.

PRODUCT USE LIMITATION

These products are intended for research use only.