

anti-PD-1 (Pembrolizumab)

Cat# A510005

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-PD-1 (Pembrolizumab)
Cat. No.	A510005
Size	1mg/5mg/20mg
Isotype	Human IgG4
Description	Pembrolizumab biosimilar is a humanized monoclonal antibody that blocks the interaction between PD-1 and its ligands, PD-L1 and PD-L2.
Species Reactivity	Human
Cas No	1374853-91-4
Host	CHO cells
Molecular Weight	149 kDa
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.
Alias	MK-3475, lambrolizumab

PRODUCT USE LIMITATION

These products are intended for research use only.