

anti-TNF-alpha (Adalimumab)

Cat# A510008

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-TNF-alpha (Adalimumab)
Cat. No.	A510008
Size	1mg/5mg/20mg
Isotype	Human IgG1
Description	Adalimumab biosimilar is a recombinant human IgG1 monoclonal antibody specific for human tumor necrosis factor (TNF), was created using phage display technology resulting in an antibody with human derived heavy and light chain variable regions and human IgG1:k constant regions.
Species Reactivity	Human
Cas No	331731-18-1
Host	CHO cells
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Molecular Weight	148 kDa
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.
Alias	LU200134, D2E7

PRODUCT USE LIMITATION

These products are intended for research use only.