

anti-CD38 (Daratumumab)

Cat# A510016

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-CD38 (Daratumumab)
Cat. No.	A510016
Size	1mg/5mg/20mg
Isotype	Human IgG1
Description	Daratumumab biosimilar is an immunoglobulin G1 kappa (IgG1κ) human monoclonal antibody against CD38 antigen, produced in a mammalian cell line (Chinese Hamster Ovary [CHO]) using recombinant DNA technology.
Species Reactivity	Human
Cas No	945721-28-8
Host	CHO cells
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Molecular Weight	148 kDa
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.
Alias	Darzalex

PRODUCT USE LIMITATION

These products are intended for research use only.