

anti-PCSK9 (Evolocumab)

Cat# A510018

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-PCSK9 (Evolocumab)
Cat. No.	A510018
Size	1mg/5mg/20mg
Isotype	Human IgG2
Description	Evolocumab biosimilar is a human monoclonal immunoglobulin G2 (IgG2) directed against human proprotein convertase subtilisin kexin 9 (PCSK9).
Species Reactivity	Human
Cas No	1256937-27-5
Host	CHO cells
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Molecular Weight	144 kDa
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.

PRODUCT USE LIMITATION

These products are intended for research use only.