

anti-EGFR (Nimotuzumab)

Cat# A510022

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-EGFR (Nimotuzumab)
Cat. No.	A510022
Size	1mg/5mg/20mg
Isotype	Chimeric IgG1
Description	Nimotuzumab biosimilar has a humanized human-mouse h-R3 heavy chain and a humanized human-mouse h-R3κ-chain. nimotuzumab binds to the epidermal growth factor receptor (EGFR), a signalling protein that normally controls cell division.
Species Reactivity	Human
Cas No	828933-51-3
Host	CHO cells
Molecular Weight	148 kDa
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.

PRODUCT USE LIMITATION

These products are intended for research use only.