

# anti-VEGFR2 (Ramucirumab)

Cat# A510029

Upon receipt, store at -80°C. Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	anti-VEGFR2 (Ramucirumab)
<b>Cat. No.</b>	A510029
<b>Size</b>	1mg/5mg/20mg
<b>Isotype</b>	human IgG1
<b>Description</b>	Ramucirumab is a recombinant human IgG1 monoclonal antibody that specifically binds to vascular endothelial growth factor receptor 2.
<b>Species Reactivity</b>	Human
<b>Cas No</b>	947687-13-0
<b>Host</b>	CHO cells
<b>Molecular Weight</b>	147 kDa
<b>Clonality</b>	Monoclonal
<b>Conjugate</b>	Non-conjugated
<b>Tested applications</b>	Functional Assay
<b>Purity</b>	>95%; < 0.75 EU/mg as determined by the LAL method
<b>Storage buffer</b>	0.01 M phosphate buffered saline (PBS) pH 7.2
<b>Storage instruction</b>	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.
<b>Alias</b>	IMC-1121B, LY3009806, IMC-1121

## PRODUCT USE LIMITATION

These products are intended for research use only.