

anti-CD20 (Rituximab)

Cat# A510030

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-CD20 (Rituximab)
Cat. No.	A510030
Size	1mg/5mg/20mg
Isotype	human IgG1
Description	Rituximab is a genetically engineered chimeric murine/human monoclonal IgG1 kappa antibody directed against the CD20 antigen.
Species Reactivity	Human
Cas No	174722-31-7
Host	CHO cells
Molecular Weight	145 kDa
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.
Alias	IDEC-C2B8

PRODUCT USE LIMITATION

These products are intended for research use only.