

anti-IL-12 and IL-23 (Ustekinumab)

Cat# A510033

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-IL-12 and IL-23 (Ustekinumab)
Cat. No.	A510033
Size	1mg/5mg/20mg
Isotype	human IgG1
Description	Ustekinumab is a human IgG1 κ monoclonal antibody against the p40 subunit of the IL-12 and IL-23 cytokines. Using DNA recombinant technology, Ustekinumab is produced in a well characterized recombinant cell line and is purified using standard bio-processing technology. The manufacturing process contains steps for the clearance of viruses.
Species Reactivity	Human
Cas No	815610-63-0
Host	CHO cells
Molecular Weight	148 kDa
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.

PRODUCT USE LIMITATION

These products are intended for research use only.