

Anti-mouse CTLA-4 (CD152)

Cat# A550001

Upon receipt, store at 4°C. Avoid repeated freeze.

INFORMATION

Product Name	Anti-mouse CTLA-4 (CD152)
Cat. No.	A550001
Size	1mg/5mg/25mg/100mg
Species Reactivity	Mouse
Immunogen	Mouse CTLA-4 IgG2a fusion protein
Concentration	3mg/ml(Size:1mg);5mg/ml(Size:5mg,25mg,100mg)
Purity	95%
Endotoxin Level	<2 EU/mg
Clone No.	UC10-4F10-11
Tested applications	in vivo/vitro CTLA-4 neutralization, FCM,WB
Molecular Weight	150 kDa
Sterility	0.2 µM filtered
Purification Method	protein A
Storage buffer	PBS buffer, pH 7.0
Storage instruction	Store the undiluted solution at 4°C in the dark to avoid freeze-thaw cycles

PRODUCT USE LIMITATION

These products are intended for research use only.