


## Anti-SARS-CoV-2 Spike Neutralizing Antibody(4A3)

Cat# A330130

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

### INFORMATION

<b>Product Name</b>	Anti-SARS-CoV-2 Spike Neutralizing Antibody(4A3)
<b>Cat. No.</b>	A330130
<b>Size</b>	50ul/100ul/200ul
<b>Immunogen</b>	Recombinant 2019-nCoV S1 Protein (mFc Tag)
<b>Host</b>	Mouse
<b>Concentration</b>	1 mg/mL
<b>Tested applications</b>	ELISA, Neutralization
<b>Application</b>	ELISA: 0.5-1 µg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification Method</b>	The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
<b>Storage buffer</b>	0.2 µm filtered solution in PBS
<b>Storage instruction</b>	Ships on ice packs. Store at -20°C
<b>Image</b>	 <p>Anti-SARS-CoV-2 Spike Neutralizing Antibody(4A3) can block Human ACE2 Protein (His Tag) and 2019-nCoV Spike Protein (RBD, mFc Tag) interaction, the IC50 for this effect is 81.5 ng/mL.</p>

### PRODUCT USE LIMITATION

These products are intended for research use only.