

Normal Rabbit Serum

Cat# P0001RbS

Store at -20°C. Avoid repeated freeze/thaw cycles.

INFORMATION

Product Name	Normal Rabbit Serum
Cat NO.	P0001Rb
Host	Rabbit
Description	Normal Serum is non-immune serum samples from rabbit that contain the usual complement of serum proteins, including the various classes of immunoglobulins that exist in healthy animals of the particular species.
Size	10 / 20 ml
Species Reactivity	Not Applicable
Tested Application	Blocking Assay, Control, ICC, IF, IHC
Application	Optimal working dilution for a specific application should be determined by the investigator.
Class	Control
Conjugate	Unconjugated
Storage Buffer	Whole Serum
Form	Liquid
Purity	Whole antiserum

PRODUCT USE LIMITATION

These products are intended for research use only.