

Anti-E (SARS-CoV) Rabbit polyclonal antibody

Cat# A140111

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Anti-E (SARS-CoV) Rabbit polyclonal antibody
Cat. No.	A140111
Size	100ug
Immunogen	Evelope protein(E)
Host	Rabbit
Clonality	Polyclonal
Tested applications	WB,ELISA,IP, IF
Conjugation	Unconjugated
Storage instruction	Store at -20°C. Avoid freeze / thaw cycle.

PRODUCT USE LIMITATION

These products are intended for research use only.