


Lamin B1 Monoclonal Antibody(7C11)

Cat# A20401PI

Upon receipt, store at -20°C . Avoid repeated freeze.

INFORMATION

Product Name	Lamin B1 Monoclonal Antibody(7C11)	
Cat. No.	A20401PI	
Size	50µl/100µl	
Uniprot	P20700	
Product type	Primary antibody	
Species Reactivity	Human,Rat,Mouse	
Immunogen	Recombinant Protein of Lamin-B1	
Host	Mouse	
Concentration	1 mg/ml	
Clonality	Monoclonal	
Tested applications	WB,IHC-p,IF,IP	
Application	WB:1/2000-5000; IP:1/200;Optimal dilutions should be determined by the end user	
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.	
Storage buffer	PBS,pH 7.4,containing 0.02% sodium azide as Preservative and 50% Glycerol	
Storage instruction	-20°C.Do not aliquot the antibody,For research use only	
Image		<p>Immunohistochemical analysis of paraffin-embedded Human-uterus tissue.</p> <ol style="list-style-type: none"> 1,Lamin B1 Monoclonal Antibody(7C11) was diluted at 1:200(4°C,overnight). 2, Sodium citrate pH 6.0 was used for antibody retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200(room temperature, 30min). Negative control was used by secondary antibody only.

PRODUCT USE LIMITATION

These products are intended for research use only.