


# GAPDH Polyclonal Antibody

Cat# A20574PI

Upon receipt, store at -20°C. Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	GAPDH Polyclonal Antibody	
<b>Cat. No.</b>	A20574PI	
<b>Size</b>	100µg/50µg	
<b>Product type</b>	Primary antibody	
<b>Species Reactivity</b>	Human,Mouse,Rat,Rb,Ch,Mk,sheep,Xenopus	
<b>Host</b>	Rabbit	
<b>Concentration</b>	1 mg/ml	
<b>Clonality</b>	Polyclonal	
<b>Tested applications</b>	WB,IHC-p	
<b>Application</b>	WB:1/3000;IHC:1/200;Optimal dilutions should be determined by the end user	
<b>Molecular Weight</b>	37KD	
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.	
<b>Storage buffer</b>	PBS,pH 7.4,containing 0.02% sodium azide as Preservative and 50% Glycerol	
<b>Storage instruction</b>	-20°C.Do not aliquot the antibody,For research use only	
<b>Image</b>		Western blot analysis of 293T (1), Rat brain (2), NIH 3T3 (3), Sheep Muscle (4), Rabbit testis (5), diluted at 1:20000.. Secondary antibody was diluted at 1:20000

## PRODUCT USE LIMITATION

These products are intended for research use only.