


# SOD2 Polyclonal Antibody

Cat# A21089PI

Storage at -20°C for 1 year

## INFORMATION

<b>Product Name</b>	SOD2 Polyclonal Antibody	
<b>Cat. No.</b>	A21089PI	
<b>Size</b>	100µg/50µg	
<b>Uniprot</b>	P04179	
<b>Product type</b>	Primary Antibody	
<b>Species Reactivity</b>	Human,Mouse,Rat	
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of human SOD-2.	
<b>Host</b>	Rb	
<b>Application</b>	WB, IHC	
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen	
<b>Storage instruction</b>	-20°C/1 year	
<b>Alias</b>	SOD2; Superoxide dismutase [Mn], mitochondrial	
<b>Image</b>		

## PRODUCT USE LIMITATION

These products are intended for research use only.