

Mox1 Polyclonal Antibody

Cat# A21372PI

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Mox1 Polyclonal Antibody
Cat. No.	A21372PI
Size	100µg/50µg/20ug
Uniprot	Human Q9Y5S8, Mouse Q8CIZ9
Product type	Primary antibody
Species Reactivity	Huma, Mouse
Immunogen	The antiserum was produced against synthesized peptide derived from human NOX1. AA range:436-485
Host	Rabbit
Concentration	1 mg/ml
Clonality	Polyclonal
Tested applications	WB,ELISA
Application	Western Blot: 1/500 - 1/2000. ELISA: 1/5000.
Molecular Weight	65 KD
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage instruction	-20°C/1 year
Alias	NOX1; MOX1; NOH1; NADPH oxidase 1; NOX-1; Mitogenic oxidase 1; MOX-1; NADH/NADPH mitogenic oxidase subunit P65-MOX; NOH-1

PRODUCT USE LIMITATION

These products are intended for research use only.