


Axl Polyclonal Antibody

Cat# A21937PI

-20°C/1 year

INFORMATION

Product Name	Axl Polyclonal Antibody
Cat. No.	A21937PI
Size	100µg/50µg/20µg
Uniprot	P30530
Product type	Primary Antibody
Species Reactivity	Human,Mouse,Rat
Immunogen	Synthesized peptide derived from Axl . at AA range: 630-710
Host	Rabbit
Concentration	1 mg/ml
Application	WB,ELISA
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage instruction	-20°C/1 year
Alias	AXL; UFO; Tyrosine-protein kinase receptor UFO; AXL oncogene
Image	 <p>293T</p> <p>170— 130— 100— Axl 70— 55— 40— 35— 25— 15—</p>

PRODUCT USE LIMITATION

These products are intended for research use only.