



IGFBP1 Polyclonal Antibody

Cat# A22477PI

Upon receipt, store at -20°C . Avoid repeated freeze.

INFORMATION

Product Name	IGFBP1 Polyclonal Antibody	
Cat. No.	A22477PI	
Size	100ug/50ug/20ug	
Uniprot	P08833	
Product type	Primary antibody	
Species Reactivity	Human	
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IGFBP1. AA range:91-140	
Host	Rabbit	
Concentration	1 mg/ml	
Tested applications	WB,IHC-p,IF(paraffin section),ELISA	
Application	IF(paraffin section) 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000.	
Molecular Weight	28KD	
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.	
Storage instruction	-20°C/1 year	
Alias	IGFBP1; IBP1; Insulin-like growth factor-binding protein 1; IBP-1; IGF-binding protein 1; IGFBP-1; Placental protein 12; PP12	
Image		Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100
		Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at



PRODUCT USE LIMITATION

These products are intended for research use only.