


IGFBP1 Polyclonal Antibody

Cat# A22842PI

Upon receipt, store at -20°C . Avoid repeated freeze.

INFORMATION

Product Name	IGFBP1 Polyclonal Antibody
Cat. No.	A22842PI
Size	100ug/50ug/20ug
Uniprot	P08833
Product type	Primary antibody
Species Reactivity	Human
Immunogen	The antiserum was produced against synthesized peptide derived from human IBP1. AA range:170-210
Host	Rabbit
Concentration	1 mg/ml
Tested applications	WB,IHC-p,IF(paraffin section),ELISA
Application	IF(paraffin section) 1:50-200 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage instruction	-20°C/1 year
Alias	IGFBP1; IBP1; Insulin-like growth factor-binding protein 1; IBP-1; IGF-binding protein 1; IGFBP-1; Placental protein 12; PP12
Image	 <p>Western Blot analysis of HepG2 cells using IGFBP1 Polyclonal Antibody</p>

PRODUCT USE LIMITATION

These products are intended for research use only.